

# GENESIS™

## TUBE SEALERS



**RAPID SEAL™  
SE850RS**  
MULTI-HEAD SEALER



**RAPID SEAL™  
SE450RS**  
MULTIPLE STATION UNIT



**RAPID SEAL™  
SE250RS**  
STANDARD BENCHTOP UNIT



**RAPID SEAL™  
SE700RS**  
PORTABLE UNIT

## Genesis Sealers... Always cool, never overheat!

- Compatible with all types of tubing
- No warm-up required – Instant On
- Notched Seal™ for easy separation
- “Intelligent” RF Board
- Lightweight, compact and durable
- Bench top, portable, multi-head, industrial devices available
- 1-year warranty parts and labor

**NEW**



**GENESIS BPS™**

Technology is in our Blood.™

# GENESIS™ TUBE SEALERS

WITH RAPID SEAL™ TECHNOLOGY



MODEL #	RAPID SEAL™ SE250RS STANDARD BENCHTOP UNIT	RAPID SEAL™ SE450RS MULTIPLE STATION UNIT	RAPID SEAL™ SE700RS PORTABLE UNIT	BIOTECH SEALER SE800 LARGER TUBING & INDUSTRIAL USE	RAPID SEAL™ SE850RS MULTI-HEAD SEALER	
<b>Workload</b>	High volume	Medium volume	Light volume/portable	Medium volume	High volume	
<b>Sealing Head Configuration</b>	Fixed sealing head	Hand-held sealing unit	Hand-held sealing unit	Heavy-duty, hand-held sealing gun	Multi-head configuration	
<b>Power Source</b>	110/220 V A/C 50/60 Hz	110/220 V A/C 50/60 Hz	12V/3 AH NiMH Battery (rechargeable)	110/220 V A/C	110/220 V A/C 50/60 Hz	
<b>Tubing Diameter</b>	2 to 6mm (1/4")	2 to 6mm (1/4")	2 to 6mm (1/4")	5 to 16mm	2 to 6mm (1/4")	
<b>DIMENSIONS</b>	<b>Main Unit</b> Inches mm	6.9W x 11.8D x 5.9H 175W x 300D x 150H	4.2W x 10.2D x 5.9H 106W x 260D x 150H	5.8W x 8.9D x 1.4H 148W x 226D x 34H	7.8W x 11.4L x 6.4H 195W x 285L x 160H	22.1W x 15D x 9.6H 561W x 380D x 244H
	<b>Remote Sealing Unit/Sealing Gun</b> Inches mm	N/A	1.1W x 8.1D x 1.8H 28W x 204D x 44H	1.1W x 8.1D x 1.8H 28W x 204D x 44H	5.4W x 1L x 5.5H 135W x 26L x 137H	N/A
<b>WEIGHT</b>	<b>Unit</b> lbs Kg	13.3 lbs 6.0 Kg	7.7 lbs 3.5 Kg	3.3 lbs 1.5 Kg	10.7 lbs 4.8 Kg	33.6 lbs 15.2 Kg
	<b>Sealing Gun</b> lbs Kg	N/A	0.5 lbs 0.2 Kg	0.5 lbs 0.2 Kg	1.1 lbs 0.5 Kg	N/A
<b>Special Features:</b>	<ul style="list-style-type: none"> <li>• RAPID SEAL™ seals in 0.8 to 1.2 seconds</li> <li>• Retractable splash guard</li> <li>• Featherlite sealing action™</li> <li>• No warm up required</li> <li>• “Intelligent” RF Board</li> </ul>	<ul style="list-style-type: none"> <li>• RAPID SEAL™ seals in 0.8 to 1.2 seconds</li> <li>• 6' standard cable</li> <li>• Easy to clean sealing head</li> <li>• No warm up required</li> </ul>	<ul style="list-style-type: none"> <li>• RAPID SEAL™ seals in 0.8 to 1.2 seconds</li> <li>• 6' standard cable</li> <li>• Lightweight canvas carry bag included</li> <li>• 6 hour fast charge battery and charger</li> <li>• No warm up required</li> </ul>	<ul style="list-style-type: none"> <li>• Biotechnology, Pharmaceutical and Industrial use</li> <li>• Adjustable sealing time</li> <li>• 6' standard cable</li> <li>• Seals tubing up to 16mm</li> <li>• No warm up required</li> </ul>	<ul style="list-style-type: none"> <li>• RAPID SEAL™ seals in 1 to 3 seconds</li> <li>• 3 1/8" between sealing heads</li> <li>• Automatic or foot pedal control</li> <li>• On/Off selectability for each head</li> <li>• Left to Right or Right to Left sealing</li> <li>• 3 color LED with alarm</li> <li>• Available with 4, 5 or 6 heads</li> <li>• No warm up required</li> </ul>	

To order, call Genesis BPS at  
**1-866-712-5663**

