

CAT# DESCRIPTION / (MFG#)

**TERUFLEX CPDA-1 BLOOD BAGS WITH
STANDARD DONOR TUBING** (Terumo)

4BBSCD456A	SINGLE BLOOD BAG (BB*SCD456A) 16 ga. ultra thinwall donor needle, 450 mL XT150 plastic container. 60/cs.
4BBDCD456A	DOUBLE BLOOD BAG (BB*DCD456A) 16 ga. ultra thinwall donor needle, 450 mL XT150 primary container with one 400 mL satellite bag. 30/cs.
4BBTCD456A	TRIPLE BLOOD BAG (BB*TCD456A) 16 ga. ultra thinwall donor needle, 450 mL XT150 primary container with one 500 mL XT612 plus 5-day platelet bag and one 400 mL satellite bag. 30/cs.

**TERUFLEX CPDA-1 BLOOD BAG SYSTEM WITH
BLOOD SAMPLING ARM** (Terumo)

For use with 4XXMN2000T and 43P-1316R (p. 28)
multi-sample luer adapter and standard tube holder.

4BSCD456A4	SINGLE BLOOD BAG (BB*SCD456A4) 16 ga. ultra thinwall donor needle, 450 mL XT150 plastic container. 48/cs.
4BDCD456A4	DOUBLE BLOOD BAG (BB*DCD456A4) 16 ga. ultra thinwall donor needle, 450 mL XT150 primary container with one 400 mL satellite bag. 30/cs.
4BTCD456A4	TRIPLE BLOOD BAG (BB*TCD456A4) 16 ga. ultra thinwall donor needle, 450 mL XT150 primary container with one 500 mL XT612 plus 5-day platelet bag and one 400 mL satellite bag. 30/cs.
4BQCD456A4	QUADRUPLE BLOOD BAG (BB*QCD456A4) 16 ga. ultra thinwall donor needle, 450 mL XT150 primary container with one 500 mL XT612 plus 5-day platelet bag and two 400 mL satellite bags. 24/cs.

CAT# DESCRIPTION / (MFG#)

**TERUFLEX CPD/OPTISOL BLOOD BAGS WITH
STANDARD DONOR TUBING** (Terumo)

4BBAGD456A	DOUBLE BLOOD BAG (BB*AGD456A) 16 ga. ultra thinwall donor needle, 450 mL XT150 primary container with one 400 mL satellite bag. 24/cs.
4BBAGT456A	TRIPLE BLOOD BAG (BB*AGT456A) 16 ga. ultra thinwall donor needle, 450 mL XT150 primary container with one 500 mL XT612 plus 5-day platelet bag and one 400 mL satellite bag. 24/cs.
4BBAGQ456A	QUADRUPLE BLOOD BAG (BB*AGQ456A) 16 ga. ultra thinwall donor needle, 450 mL XT150 primary container with one 500 mL XT612 plus 5-day platelet bag and two 400 mL satellite bag. 24/cs.

**TERUFLEX CPD/OPTISOL BLOOD BAG SYSTEM WITH
BLOOD SAMPLING ARM** (Terumo)

For use with 4XXMN2000T and 43P-1316R.

4BAGD456A4	DOUBLE BLOOD BAG (BB*AGD456A4) 16 ga. ultra thinwall donor needle, 450 mL XT150 primary container with one 400 mL satellite bag. 24/cs.
4BAGT456A4	TRIPLE BLOOD BAG (BB*AGT456A4) 16 ga. ultra thinwall donor needle, 450 mL XT150 primary container with one 500 mL XT612 plus 5-day platelet bag and one 400 mL satellite bag. 24/cs.
4BAGQ456A4	QUADRUPLE BLOOD BAG (BB*AGQ456A4) 16 ga. ultra thinwall donor needle, 450 mL XT150 primary container with one 500 mL XT612 plus 5-day platelet bag and two 400 mL satellite bags. 24/cs.

CAT# DESCRIPTION / (MFG#)

**TERUMO (TERUFLEX) BLOOD BAGS WITH SAMPLING ARM®
AND DONOR CARE® NEEDLE GUARD PRE ATTACHED:**

450mL CPDA-1:

- 4BBSCD456A3 SINGLE BLOOD BAG
16 gauge ultra-thinwall needle, a 450mL XT150
primary container. 36/case
- 4BBDCD456A3 DOUBLE BLOOD BAG
16 gauge ultra-thinwall needle, a 450mL XT150
primary container, one 400mL satellite bag. 30/case
- 4BBQCD456A3 QUADRUPLE BLOOD BAG
16 gauge ultra-thinwall needle, a 450mL XT150
primary container, two 400mL satellite bags and
one XT-612 platelet storage bag. 24/case

450mL CPD/OPTISOL:

- 4BBAGD456A3 DOUBLE BLOOD BAG
16 gauge ultra-thinwall needle, a 450mL XT150
primary container, one 400mL satellite bag. 24/case
- 4BBAGT456A3 TRIPLE BLOOD BAG
16 gauge ultra-thinwall needle, a 450mL XT150
primary container, one 400mL satellite bag and
one XT-612 platelet storage bag. 24/case
- 4BBAGQ456A3 QUADRUPLE BLOOD BAG
16 gauge ultra-thinwall needle, a 450mL XT150
primary container, two 400mL satellite bags and
one XT-612 platelet storage bag.
(Special order; Minimum order required)

500mL CPDA-1:

- 4BBSCD506A3 SINGLE BLOOD BAG
16 gauge ultra-thinwall needle, a 500mL XT150
primary container. 36/case
- 4BBDCD506A3 DOUBLE BLOOD BAG
16 gauge ultra-thinwall needle, a 500mL XT150
primary container, one 400mL satellite bag. 30/case
- 4BBTCD506A3 TRIPLE BLOOD BAG
16 gauge ultra-thinwall needle, a 500mL XT150
primary container, one 400mL satellite bag and
one XT-612 platelet storage bag.
(Special order; Minimum order required)
- 4BBQCD506A3 QUADRUPLE BLOOD BAG
16 gauge ultra-thinwall needle, a 500mL XT150
primary container, two 400mL satellite bags and
one XT-612 platelet storage bag. 24/case

500mL CPD/OPTISOL:

- 4BBAGD506A3 DOUBLE BLOOD BAG
16 gauge ultra-thinwall needle, a 500mL XT150
primary container, one 400mL satellite bag. 24/case
- 4BBAGT506A3 TRIPLE BLOOD BAG
16 gauge ultra-thinwall needle, a 500mL XT150
primary container, one 400mL satellite bag and
one XT-612 platelet storage bag. 24/case

CAT# DESCRIPTION / (MFG#)

**TERUMO (TERUFLEX) BLOOD BAGS WITH DIVERSION
BLOOD SAMPLING ARM®:**

450mL CPD/OPTISOL:

- 4BBAGD456A2 DOUBLE BLOOD BAG
16 gauge ultra-thinwall needle, a 450mL XT150
primary container, one 400mL satellite bag. 24/case
- 4BBAGT456A2 TRIPLE BLOOD BAG
16 gauge ultra-thinwall needle, a 450mL XT150
primary container, one 400mL satellite bag and
one XT-612 platelet storage bag. 24/case
- 4BBAGQ456A2 QUADRUPLE BLOOD BAG
16 gauge ultra-thinwall needle, a 450mL XT150
primary container, two 400mL satellite bags and
one XT-612 platelet storage bag. 18/case

500mL CPD/OPTISOL:

- 4BBAGD506A2 DOUBLE BLOOD BAG
16 gauge ultra-thinwall needle, a 500mL XT150
primary container, one 400mL satellite bag. 24/case
- 4BBAGT506A2 TRIPLE BLOOD BAG
16 gauge ultra-thinwall needle, a 500mL XT150
primary container, one 400mL satellite bag and
one XT-612 platelet storage bag. 24/case
- 4BBAGQ506A2 QUADRUPLE BLOOD BAG
16 gauge ultra-thinwall needle, a 500mL XT150
primary container, two 400mL satellite bags and
one XT-612 platelet storage bag. 18/case

**TERUMO IMUFLEX® WB-RP BLOOD BAG WITH
IN-LINE LEUKOCYTE REDUCTION FILTER AND BLOOD
SAMPLING ARM®:**

500mL CPD/OPTISOL:

- 4BBWQ506A3 TRIPLE BLOOD BAG
contains an in-line IMUFLEX® leukocyte reduction
filter connected to a 500mL XT150 primary container,
two 400mL satellite bags. Donor Care® Needle
Guard is pre-attached. 18/case