



**GENESIS BPS™**

Technology is in our Blood.<sup>SM</sup>

# CROSS REFERENCE GUIDE

## BLOOD BAGS

| <b>CPDA-1</b>  |        | <b>Baxter</b> | <b>Baxter</b> | Ⓣ <b>TERUMO</b> | Ⓣ <b>TERUMO</b> |
|--|--------|---------------|---------------|-----------------|-----------------|
|  |        | <b>450 mL</b> | <b>500 mL</b> | <b>450 mL</b>   | <b>500 mL</b>   |
| Standard tubing  | single | –             | –             | 4BBSCD456A      | –               |
|  | double | –             | –             | 4BBDCD456A      | –               |
|  | triple | –             | –             | 4BBTCD456A      | –               |
|  | quad   | –             | –             | –               | –               |
| Blood Sampling Arm   | single | –             | –             | 4BBSCD456A4     | –               |
|  | double | –             | –             | 4BBDCD456A4     | –               |
|  | triple | –             | –             | 4BBTCD456A4     | –               |
|  | quad   | –             | –             | 4BBQCD456A4     | –               |
| Lab-site sampling arm w / ITL donor guard needle protector | single | 4R3610NM      | 4R3615NM      | 4BBSCD456A3     | 4BBSCD506A3     |
|  | double | 4R3620NM      | 4R3625NM      | 4BBDCD456A3     | 4BBDCD506A3     |
|  | triple | 4R3631NM      | 4R3635NM      | 4BBTCD456A3     | –               |
|  | triple | 4R3630NM      | 4R3636NM      | –               | –               |
|  | quad   | 4R3640NM      | 4R3645NM      | 4BBQCD456A3     | 4BBQCD506A3     |
|  | quint  | –             | 4R3655NM      | –               | –               |

| <b>CPD AND ADDITIVE SOLUTIONS</b>                          |        | <b>Adsol (CPD)</b> | <b>Adsol (CPD)</b> | <b>Optisol (CPD)</b> | <b>Optisol (CPD)</b> |
|--|--------|--------------------|--------------------|----------------------|----------------------|
|  |        | <b>450 mL</b>      | <b>500 mL</b>      | <b>450 mL</b>        | <b>500 mL</b>        |
| Standard tubing  | double | –                  | –                  | 4BBAGD456A           | –                    |
|  | triple | –                  | –                  | 4BBAGT456A           | –                    |
|  | quad   | –                  | –                  | 4BBAGQ456A           | –                    |
| Blood Sampling Arm   | double | –                  | –                  | 4BBAGD456A4          | –                    |
|  | triple | –                  | –                  | 4BBAGT456A4          | –                    |
|  | quad   | –                  | –                  | 4BBAGQ456A4          | –                    |
| Lab-site sampling arm w / ITL donor guard needle protector | single | –                  | –                  | –                    | –                    |
|  | double | 4R3420NM           | 4R3425NM           | 4BBAGD456A3          | 4BBAGD506A3          |
|  | triple | 4R3431NM           | 4R3435NM           | 4BBAGT456A3          | 4BBAGT506A3          |
|  | triple | –                  | 4R3436NM           | –                    | –                    |
|  | quad   | 4R3441NM           | 4R3446NM           | –                    | –                    |
|  | quad   | –                  | 4R3448NM           | –                    | –                    |

**Genesis BPS** carries most blood bags, transfer bags, blood I.V. sets and transfer sets from all major manufacturers. If you do not see the item you need above, please ask your representative or our customer service department for a price quote.





# CROSS REFERENCE GUIDE


## BLOOD BAGS


**GENESIS BPS™**

Technology is in our Blood.™

|  <b>TERUMO</b><br>(TERUFLEX) WITH SAMPLING ARM | CPDA-1<br>450 mL | CPD/OPTISOL<br>450 mL |
|---|------------------|-----------------------|
|   | single           | 4BBSCD456A4           |
| double  | 4BBDCD456A4      | 4BBAGD456A4           |
| triple  | 4BBTCD456A4      | 4BBAGT456A4           |
| quad  | 4BBQCD456A4      | 4BBAGQ456A4           |

|  <b>TERUMO</b><br>(TERUFLEX) BLOOD BAGS W/ DIVERSION<br>BLOOD SAMPLING ARM | CPD/OPTISOL<br>450 mL | CPD/OPTISOL<br>500 mL |
|---|-----------------------|-----------------------|
|   | double                | 4BBAGD456A2           |
| triple  | 4BBAGT456A2           | 4BBAGT506A2           |
| quad  | 4BBAGQ456A2           | 4BBAGQ506A2           |

|  <b>TERUMO</b><br>(TERUFLEX) CPD/OPTISOL 250mL/ 17G WITH SAMPLING ARM | 4BBAGQ257A4 |
|--|-------------|
|  |             |

|  <b>TERUMO</b><br>(TERUFLEX) CPD/OPTISOL DOUBLE W/ SAMPLING ARM 500mL<br>W/ IN LINE IMUFLEX FILTER AND DONOR GUARD NEEDLE PROTECTOR | 4BB*WQQ506A3 |
|--|--------------|
|  |              |

Genesis BPS carries most blood bags, transfer bags, blood I.V. sets and transfer sets from all major manufacturers. If you do not see the item you need above, please ask your representative or our customer service department for a price quote.



**GENESIS** BPS™

Technology is in our Blood.™

## CROSS REFERENCE GUIDE

### LEUKOCYTE REDUCTION BLOOD BAGS TRANSFER BAGS

| <b>Baxter</b> LEUKOCYTE REDUCTION BLOOD BAGS                   | 450 mL   | 500 mL   |
|--|----------|----------|
| CPDA-1 / RS 2000 Integral Filter - Double - Sampling Site      | 4R3310NM | 4R3315NM |
| CPD / ADSOL / RZ 2000 Integral Filter - Double - Sampling Site | 4R3316NM | 4R3317NM |
| CPD / ADSOL / RZ 2000 Integral Filter - Triple - Sampling Site | –        | 4R3319NM |
| CPD / ADSOL / RS 2000 Integral Filter - Double - Sampling Site | 4R3303NM | 4R3308NM |
| CPD / ADSOL / RS 2000 Integral Filter -Triple - Sampling Site  | –        | 4R3309NM |

| TRANSFER BAGS                                      | GENESIS BPS™ | <b>Baxter</b> | Ⓣ TERUMO   |
|--|--------------|---------------|------------|
| 75 mL Pedi-Pak Quad with filter                    | 402-04F      | –             | –          |
| 75 mL Pedi-Pak Quad without filter                 | 402-04       | –             | –          |
| 75 mL Pedi-Pak Quad needleless Y site              | 402-01       | –             | –          |
| 75 mL Pedi-Pak Single                              | 402-024      | –             | –          |
| 150 mL Single                                      | –            | 4R2001        | 4BT015CB70 |
| 150 mL Quadruple                                   | –            | 4R2004        | –          |
| 300 mL Single with coupler                         | –            | 4R2014        | 4BT030CB71 |
| 600 mL Single with coupler                         | –            | 4R2023        | 4BT060CB71 |
| 600 mL Single with needle adaptor                  | –            | 4R2024        | –          |
| 800 mL Single with 2 couplers without outlet ports | –            | 4R2054        | –          |
| 800 mL Single with 2 couplers & outlet ports       | –            | 4R2055        | 4BT080BB71 |
| 1000 mL Single with coupler                        | –            | 4R2032        | 4BT100BB71 |
| 2000 mL Single with coupler                        | –            | 4R2041        | 4BT200BB71 |
| Platelet Pooling Unit 600 mL / 8 couplers          | –            | 4R2027        | –          |
| Single dry bag w/16g needle 600 mL                 | –            | –             | 4BBD606A   |
| Single dry bag w/15g needle 600 mL                 | –            | 4R0001        | –          |

**Genesis BPS** carries Leukocyte reduction products.  
Please call our customer service department for a price quote.

| CAT#   | DESCRIPTION   |
|--|---|
| <b>BLOOD-PACK CONTAINER SYSTEMS WITH CPDA-1 (Baxter)</b> |   |
|  | <b>Lab-Site Sampling Site (use with SampLink™ Access Device 4R4215) ITL Donorguard needle protector.</b><br>450 mL primary pack size.   |
| 4R3610NM   | SINGLE BLOOD-PACK UNIT<br>450 mL PL 146 plastic Blood-Pack unit, 16 ga. needle. 60/case   |
| 4R3620NM   | DOUBLE BLOOD-PACK UNIT<br>450 mL PL 146 plastic primary container, 16 ga. needle,<br>with one empty 400 mL PL 146 plastic Transfer Pack. 48/case.   |
| 4R3630NM   | TRIPLE BLOOD-PACK UNIT<br>450 mL PL 146 plastic primary container, 16 ga. needle,<br>with two empty 400 mL PL 146 plastic Transfer Packs. 36/case.  |
| 4R3631NM   | TRIPLE BLOOD-PACK UNIT<br>450 mL PL 146 plastic primary container, 16 ga. needle,<br>with one empty 300 mL PL 732 plastic Transfer Pack and<br>one 400 PL 146 plastic Transfer Pack. 36/case.     |
| 4R3640NM   | QUADRUPLE BLOOD-PACK UNIT<br>450 mL PL 146 plastic primary container, 16 ga. needle,<br>with three empty 400 mL PL 146 plastic Transfer Packs. 24/case.   |
|  | <b>Lab-Site Sampling Site (use with SampLink™ Access Device 4R4215) and ITL Donorguard needle protector.</b><br>500 mL primary container.   |
| 4R3615NM   | SINGLE BLOOD-PACK UNIT<br>500 mL PL 146 plastic Blood-Pack unit, 16 ga. needle. 60/case.  |
| 4R3625NM   | DOUBLE BLOOD-PACK UNIT<br>500 mL PL 146 plastic primary container, 16 ga. needle,<br>with one empty 400 mL PL 146 plastic Transfer Pack. 48/case.   |
| 4R3635NM   | TRIPLE BLOOD-PACK UNIT<br>500 mL PL 146 plastic primary container, 16 ga. needle,<br>with two empty 400 mL PL 146 plastic Transfer Packs. 36/case.  |
| 4R3636NM   | TRIPLE BLOOD-PACK UNIT<br>500 mL PL 146 plastic primary container, 16 ga. needle,<br>with one empty 300 mL PL 732 plastic Transfer Pack, and<br>one 400 mL PL 146 plastic Transfer Pack. 36/case. |
| 4R3645NM   | QUADRUPLE BLOOD-PACK UNIT<br>500 mL PL 146 plastic primary container, 16 ga. needle, with<br>three empty 400 mL PL 146 plastic Transfer Packs. 24/case.   |
| 4R3655NM   | QUINTUPLE BLOOD-PACK UNIT<br>500 mL PL 146 plastic primary container, 15 ga. needle,<br>with four empty 400 mL PL 146 plastic Transfer Packs. 24/case.  |

| CAT#  | DESCRIPTION  |
|---|--|
| <b>BLOOD-PACK CONT. SYST. WITH CPD/ADSOL (Baxter)</b> |  |
|   | <b>Lab-Site Sampling Site (use with SampLink™ Access Device 4R4215) and ITL Donorguard needle protector.</b><br>450 mL primary pack size.  |
| 4R3420NM  | DOUBLE BLOOD-PACK UNIT<br>450 mL PL 146 plastic primary container, 16 ga. needle,<br>with one empty 400 mL PL 146 plastic Transfer Pack<br>and one 400 mL PL 146 plastic satellite pack<br>containing 100 mL Adsol solution. 30/case   |
| 4R3431NM  | TRIPLE BLOOD-PACK UNIT<br>450 mL PL 146 plastic primary container, 16 ga. needle,<br>with one empty 300 mL PL 732 plastic Transfer Pack,<br>and one 400 mL PL 146 plastic Transfer Pack with<br>100 mL Adsol solution. 30/case.  |
| 4R3441NM  | QUADRUPLE BLOOD-PACK UNIT<br>450 mL PL 146 plastic primary container, 16 ga. needle,<br>with one empty 300 mL PL 732 plastic Transfer Pack,<br>one empty 400 mL PL 146 plastic Transfer Pack, and<br>one 400 mL PL 146 plastic Transfer Pack with 100 mL<br>Adsol solution. 24/case.   |
|   | <b>Lab-Site Sampling Site (use with SampLink™ Access Device 4R4215) and ITL Donorguard needle protector.</b><br>(500 mL primary container.)  |
| 4R3425NM  | DOUBLE BLOOD-PACK UNIT<br>500 mL PL 146 plastic primary container, 16 ga. needle,<br>with one empty 400 mL PL 146 plastic Transfer Pack and<br>one 400 mL PL 146 plastic nonreusable satellite pack<br>containing 110 mL Adsol solution. 30/case.  |
| 4R3435NM  | TRIPLE BLOOD-PACK UNIT<br>500 mL PL 146 plastic primary container, 16 ga. needle,<br>with one empty 400 mL PL 146 plastic Transfer Pack, and<br>one 400 mL PL 146 plastic Transfer Pack with 110 mL Adsol<br>solution. 30/case.  |
| 4R3436NM  | TRIPLE BLOOD-PACK UNIT<br>500 mL PL 146 plastic primary container, 16 ga. needle,<br>with one empty 300 mL PL 732 plastic Transfer Pack, and<br>one 400 mL PL 146 plastic Transfer Pack with 110 mL Adsol<br>solution. 30/case.  |
| 4R3446NM  | QUADRUPLE BLOOD-PACK UNIT<br>w/Labsite Sampling Site, Donor Care Needle Guard,<br>16ga needle, 500mL PL 146 plastic primary container,<br>one empty 300mL PL 732 plastic Transfer Pack unit,<br>one empty 400mL PL1240 plastic Transfer Pack unit, and<br>one 400mL PL 146 plastic Transfer Pack unit 2/110mL<br>ADSOL solution. 24/case |
| 4R3448NM  | QUADRUPLE BLOOD-PACK UNIT<br>w/Labsite Sampling Site, Donor Care Needle Guard,<br>16ga needle, 500mL PL 146 plastic primary container,<br>two empty 400mL PL 1240 plastic Transfer Pack units,<br>and one 400mL PL 146 plastic Transfer Pack unit w/110mL<br>ADSOL Solution. 24/case   |
| 4R4215  | SAMPLINK™ ACCESS DEVICE<br>for use with Blood-Pack units with Lab-Site Sampling Site.<br>60/case.  |

**CAT# DESCRIPTION / (MFG#)**

**TERUFLEX CPDA-1 BLOOD BAGS WITH  
STANDARD DONOR TUBING** (Terumo)

|            |  |
|------------|--|
| 4BBSCD456A | SINGLE BLOOD BAG (BB*SCD456A)<br>16 ga. ultra thinwall donor needle, 450 mL XT150 plastic container. 60/cs.  |
| 4BBDCD456A | DOUBLE BLOOD BAG (BB*DCD456A)<br>16 ga. ultra thinwall donor needle, 450 mL XT150 primary container with one 400 mL satellite bag. 30/cs.  |
| 4BBTCD456A | TRIPLE BLOOD BAG (BB*TCD456A)<br>16 ga. ultra thinwall donor needle, 450 mL XT150 primary container with one 500 mL XT612 plus 5-day platelet bag and one 400 mL satellite bag. 30/cs. |

**TERUFLEX CPDA-1 BLOOD BAG SYSTEM WITH  
BLOOD SAMPLING ARM** (Terumo)

For use with 4XXMN2000T and 43P-1316R (p. 28)  
multi-sample luer adapter and standard tube holder.

|            |   |
|------------|---|
| 4BSCD456A4 | SINGLE BLOOD BAG (BB*SCD456A4)<br>16 ga. ultra thinwall donor needle, 450 mL XT150 plastic container. 48/cs.  |
| 4BDCD456A4 | DOUBLE BLOOD BAG (BB*DCD456A4)<br>16 ga. ultra thinwall donor needle, 450 mL XT150 primary container with one 400 mL satellite bag. 30/cs.  |
| 4BTCD456A4 | TRIPLE BLOOD BAG (BB*TCD456A4)<br>16 ga. ultra thinwall donor needle, 450 mL XT150 primary container with one 500 mL XT612 plus 5-day platelet bag and one 400 mL satellite bag. 30/cs.     |
| 4BQCD456A4 | QUADRUPLE BLOOD BAG (BB*QCD456A4)<br>16 ga. ultra thinwall donor needle, 450 mL XT150 primary container with one 500 mL XT612 plus 5-day platelet bag and two 400 mL satellite bags. 24/cs. |

**CAT# DESCRIPTION / (MFG#)**

**TERUFLEX CPD/OPTISOL BLOOD BAGS WITH  
STANDARD DONOR TUBING** (Terumo)

|            |   |
|------------|---|
| 4BBAGD456A | DOUBLE BLOOD BAG (BB*AGD456A)<br>16 ga. ultra thinwall donor needle, 450 mL XT150 primary container with one 400 mL satellite bag. 24/cs.   |
| 4BBAGT456A | TRIPLE BLOOD BAG (BB*AGT456A)<br>16 ga. ultra thinwall donor needle, 450 mL XT150 primary container with one 500 mL XT612 plus 5-day platelet bag and one 400 mL satellite bag. 24/cs.    |
| 4BBAGQ456A | QUADRUPLE BLOOD BAG (BB*AGQ456A)<br>16 ga. ultra thinwall donor needle, 450 mL XT150 primary container with one 500 mL XT612 plus 5-day platelet bag and two 400 mL satellite bag. 24/cs. |

**TERUFLEX CPD/OPTISOL BLOOD BAG SYSTEM WITH  
BLOOD SAMPLING ARM** (Terumo)

For use with 4XXMN2000T and 43P-1316R.

|            |   |
|------------|---|
| 4BAGD456A4 | DOUBLE BLOOD BAG (BB*AGD456A4)<br>16 ga. ultra thinwall donor needle, 450 mL XT150 primary container with one 400 mL satellite bag. 24/cs.  |
| 4BAGT456A4 | TRIPLE BLOOD BAG (BB*AGT456A4)<br>16 ga. ultra thinwall donor needle, 450 mL XT150 primary container with one 500 mL XT612 plus 5-day platelet bag and one 400 mL satellite bag. 24/cs.     |
| 4BAGQ456A4 | QUADRUPLE BLOOD BAG (BB*AGQ456A4)<br>16 ga. ultra thinwall donor needle, 450 mL XT150 primary container with one 500 mL XT612 plus 5-day platelet bag and two 400 mL satellite bags. 24/cs. |

CAT# DESCRIPTION / (MFG#)

**TERUMO (TERUFLEX) BLOOD BAGS WITH SAMPLING ARM®  
AND DONOR CARE® NEEDLE GUARD PRE ATTACHED:**

**450mL CPDA-1:**

|             |   |
|-------------|---|
| 4BBSCD456A3 | SINGLE BLOOD BAG<br>16 gauge ultra-thinwall needle, a 450mL XT150 primary container. 36/case  |
| 4BBDCD456A3 | DOUBLE BLOOD BAG<br>16 gauge ultra-thinwall needle, a 450mL XT150 primary container, one 400mL satellite bag. 30/case   |
| 4BBQCD456A3 | QUADRUPLE BLOOD BAG<br>16 gauge ultra-thinwall needle, a 450mL XT150 primary container, two 400mL satellite bags and one XT-612 platelet storage bag. 24/case |

**450mL CPD/OPTISOL:**

|             |  |
|-------------|--|
| 4BBAGD456A3 | DOUBLE BLOOD BAG<br>16 gauge ultra-thinwall needle, a 450mL XT150 primary container, one 400mL satellite bag. 24/case  |
| 4BBAGT456A3 | TRIPLE BLOOD BAG<br>16 gauge ultra-thinwall needle, a 450mL XT150 primary container, one 400mL satellite bag and one XT-612 platelet storage bag. 24/case  |
| 4BBAGQ456A3 | QUADRUPLE BLOOD BAG<br>16 gauge ultra-thinwall needle, a 450mL XT150 primary container, two 400mL satellite bags and one XT-612 platelet storage bag.<br>(Special order; Minimum order required) |

**500mL CPDA-1:**

|             |  |
|-------------|--|
| 4BBSCD506A3 | SINGLE BLOOD BAG<br>16 gauge ultra-thinwall needle, a 500mL XT150 primary container. 36/case   |
| 4BBDCD506A3 | DOUBLE BLOOD BAG<br>16 gauge ultra-thinwall needle, a 500mL XT150 primary container, one 400mL satellite bag. 30/case  |
| 4BBTCD506A3 | TRIPLE BLOOD BAG<br>16 gauge ultra-thinwall needle, a 500mL XT150 primary container, one 400mL satellite bag and one XT-612 platelet storage bag.<br>(Special order; Minimum order required) |
| 4BBQCD506A3 | QUADRUPLE BLOOD BAG<br>16 gauge ultra-thinwall needle, a 500mL XT150 primary container, two 400mL satellite bags and one XT-612 platelet storage bag. 24/case                                |

**500mL CPD/OPTISOL:**

|             |   |
|-------------|---|
| 4BBAGD506A3 | DOUBLE BLOOD BAG<br>16 gauge ultra-thinwall needle, a 500mL XT150 primary container, one 400mL satellite bag. 24/case                                     |
| 4BBAGT506A3 | TRIPLE BLOOD BAG<br>16 gauge ultra-thinwall needle, a 500mL XT150 primary container, one 400mL satellite bag and one XT-612 platelet storage bag. 24/case |

CAT# DESCRIPTION / (MFG#)

**TERUMO ( TERUFLEX ) BLOOD BAGS WITH DIVERSION  
BLOOD SAMPLING ARM®:**

**450mL CPD/OPTISOL:**

|             |   |
|-------------|---|
| 4BBAGD456A2 | DOUBLE BLOOD BAG<br>16 gauge ultra-thinwall needle, a 450mL XT150 primary container, one 400mL satellite bag. 24/case   |
| 4BBAGT456A2 | TRIPLE BLOOD BAG<br>16 gauge ultra-thinwall needle, a 450mL XT150 primary container, one 400mL satellite bag and one XT-612 platelet storage bag. 24/case     |
| 4BBAGQ456A2 | QUADRUPLE BLOOD BAG<br>16 gauge ultra-thinwall needle, a 450mL XT150 primary container, two 400mL satellite bags and one XT-612 platelet storage bag. 18/case |

**500mL CPD/OPTISOL:**

|             |   |
|-------------|---|
| 4BBAGD506A2 | DOUBLE BLOOD BAG<br>16 gauge ultra-thinwall needle, a 500mL XT150 primary container, one 400mL satellite bag. 24/case   |
| 4BBAGT506A2 | TRIPLE BLOOD BAG<br>16 gauge ultra-thinwall needle, a 500mL XT150 primary container, one 400mL satellite bag and one XT-612 platelet storage bag. 24/case     |
| 4BBAGQ506A2 | QUADRUPLE BLOOD BAG<br>16 gauge ultra-thinwall needle, a 500mL XT150 primary container, two 400mL satellite bags and one XT-612 platelet storage bag. 18/case |

**TERUMO IMUFLEX® WB-RP BLOOD BAG WITH  
IN-LINE LEUKOCYTE REDUCTION FILTER AND BLOOD  
SAMPLING ARM®:**

**500mL CPD/OPTISOL:**

|            |   |
|------------|---|
| 4BBWQ506A3 | TRIPLE BLOOD BAG<br>contains an in-line IMUFLEX® leukocyte reduction filter connected to a 500mL XT150 primary container, two 400mL satellite bags. Donor Care® Needle Guard is pre-attached. 18/case |
|------------|---|