

# GENESIS™

## BLOOD COLLECTION MIXER

MODEL# CM735A  
CAT# 425-CM735A



**Our mixer stands alone in the field of blood collection technology – the most accurate, durable and reliable automated mixer.**

- Weight accuracy +/- 1 gram
- Monitors blood flow
- Automatic taring
- Reduces QNS rate and over bleeds
- Rechargeable battery
- Perfect for mobile and fixed sites
- Enhanced digital display



**GENESIS BPS™**

Technology is in our Blood.™

WHEN YOU'RE OUT FOR BLOOD...TAKE US ALONG.

# GENESIS™ BLOOD COLLECTION MIXER

MODEL# CM735A

Perfect for fixed and mobile collections  
Rechargeable battery or A/C powered



## SPECIAL FEATURES

- Rocking action mixes blood with anticoagulant
- Automatically clamps tubing when collected volume reaches pre-set volume
- Enhanced Digital Displays:
  - Pre-set volume by operator
  - Collected volume in mL
  - Total weight in grams
  - Elapsed time of donation
- Audible/Visual Alarms:
  - Slow bleed
  - Low battery
  - End of collection
- Tray:
  - Deeper – Holds all bags & In-line filters
  - Magnetic base ensures stability
  - Leak-proof
  - Tubing guides
- Reduces over and under bleeds (QNS units)
- Automatically tares collection bag & filters
- Rechargeable Ni-Cad battery or 110V/220V A/C power
- Individual recharger included
- Designed for mobile use/transportation case included
- Effective for use in collecting cord blood
- UL/CSA approved
- One year limited warranty

## MODEL# CM735A SPECIFICATIONS

External Dimensions ..... 11.5"L x 8.0"W x 5.5"H  
Collection Volume ..... 0–1199 mL  
Weight Display ..... 0–1199 grams  
Elapsed Time Display ..... 0–19 minutes 59 seconds  
Target Volume Range ..... 0–1199 mL Default Value 450 mL  
or 500 mL with 5 mL increments

### Accuracy:

Weight Accuracy ..... +/- 1 gram  
Weight Tolerance ..... 0 grams to + 3 grams

### LCD Custom Designed Display

Power Requirements ..... 110V–120V A/C, 220V–240V A/C,  
50–60 Hz

Battery Pack ..... Rechargeable Ni-Cad 12V DC

### Weight:

Mixer with battery ..... 7 lbs  
Mixer with case and all accessories ..... 12 lbs  
Shipping Weight ..... 17 lbs  
Packing ..... One unit per carton

### Accessories Included:

Ni-Cad Battery Pack  
Battery Recharger Unit  
A/C Power Adapter\*  
User Manual  
Carrying/Transport Case  
Mixer Stand

### Optional Accessories:

Multi-Battery Charger  
Extended Service Agreement

\*Suitable for 110V–120V/220V–240V A/C dependent on power supply.



## Genesis™ Multi-Charger with Rapid Charge-Discharge Feature

CAT# 425-BC162

- Charge up to 6 batteries simultaneously
- 1 slot for rapid discharge/charge feature
- For maximum battery life, use the rapid discharge feature on each battery once per year.

To order, call Genesis BPS at

**1-866-712-5663**



**GENESIS BPS™**

Technology is in our Blood.™