

GENESIS ACCESSORIES



LABath Clean™

CAT# 402-05

- Water Bath Clarifier
- 240 mL concentrate



Plasma Guard Thawing Bags™

CAT# 401-03

- Protects ports during thawing process
- Unique double seal prevents/eliminates leaks



Sure Seal Thawing Bags™

CAT# 401-021

- Protects ports during thawing process
- Unique double seal prevents/eliminates leaks



4" Blood Shoe™ Blood Bag Holders

CAT# 400-354-A

- Positions blood bags upright
- Labels are legible
- Clear, durable and easily cleaned



6" Blood Shoe™ Blood Bag Holders

CAT# 400-356-A

- Positions blood bags upright
- Labels are legible
- Clear, durable and easily cleaned



GENESIS Tube Stripper

CAT# 4A-HS-002

- High Quality aluminum alloy
- Light weight and sturdy
- Ergonomically designed

Ancillary Blood Collection Supplies

GenesisBPS carries a full line of equipment and supplies to meet all your collection needs, including:

- Blood collection bags/transfer bags
- Sealing clips/grommets
- Hand/phlebotomy squeezers

GENESIS BPS™

GenesisBPS meets client needs worldwide, for the collection, processing, storage, and preparation of blood and other biological fluids.

We do so by combining advanced technology, quality standards and operational efficiencies, working together to provide products and services of outstanding precision and reliability.

GenesisBPS is positioned at the forefront of what we can do and at the leading edge of what we will do, providing ideas, solutions, and innovation to the world.

www.GenesisBPS.com

1-866-712-5663



GENESIS™

BLOOD PRODUCT SPECIALTIES

A Unique Line of Supplies for
the Blood Bank & Laboratory



**Blood Bank Specialty
Equipment and Supplies**

GENESIS BPS™

Technology is in our Blood.™

www.GenesisBPS.com
1-866-712-5663

COLLECTION

Improved Design & Features
• Battery or A/C Powered



GENESIS™ BLOOD COLLECTION MIXER CAT# 425-CM735A

This compact mixer can be operated by A/C current or rechargeable battery pack. It features electronic monitoring with digital display and conveniently automates mixing, weighing and blood flow.

- Weight accuracy +/- 1 gram
- Monitors blood flow
- Automatic taring
- Reduces QNS rate and over bleeds
- Rechargeable battery
- Perfect for mobile and fixed sites
- Enhanced digital display

GENESIS™ MULTI-CHARGER with Rapid Charge/Discharge Feature CAT# 425-BC162

- Charge up to 6 batteries simultaneously
- 1 slot for fast charge/discharge feature
- For maximum battery life, use the rapid discharge feature on each battery once per year.



GENESIS™ CM350 Blood Collection Monitor CAT# 425-CM350

- Electronically monitors and measures blood flow and volume
- Automatically clamps with pre-set volume
- Automatic weight calibration • Dual function keypad
- Rechargeable NiMH cell 7.2V/1.8ah battery



IV ADMIN SETS

- CAT# 60150-A (Vented)
- CAT# 60151C-A (Pediatric)
- CAT# 60160-A (Non-Vented)
- CAT# 60150C-A (Universal)
- CAT# 60153C-A (Universal with 15 micron filter)

- Non-DEHP
- Leur Slip and Roller Clamps

STORAGE



TCD® STERILE TUBE WELDER CAT# 480-B40-00-1245

- Small footprint
- Single wafer loading—no jamming
- Bent tube technology
- Greater wafer mass
- No validation test kits required

PEDI-PAK® PEDIATRIC 75 mL TRANSFER PACK QUAD OR SINGLE

- CAT# 402-04 (QUAD without Filter)
- CAT# 402-04F (QUAD with Filter)
- CAT# 402-24 (QUAD Needleless "Y" Site)
- CAT# 402-01 (SINGLE)

PEDI-PAK® 75 mL Transfer Packs are the ideal solution for splitting an adult size unit of blood into four equal pediatric size units.

- Maximize usage of adult size blood units
- Reduces risk of donor exposure to neonates
- Ideal for aliquotting units of rare blood and FFP
- Compatible with sterile connecting technology



PEDI-PAK® PEDI-SYRINGE FILTER™ 30 mL & 60 mL

- CAT# 405-535MJ (35 mL)
- CAT# 405-560MJ (60 mL)
- CAT# 405-530BD (30 mL)
- CAT# 405-560BD (60 mL)

- 30 mL or 60 mL sets available with 150 micron filter
- Limits neonate exposure to multiple donors
- Maximize use of blood supply
- Compatible with sterile connecting technology
- Needleless technology
- Compatible with all syringe pumps
- DEHP Free



PROCESSING



GENESIS™ TUBE SEALERS

Innovative multi-head designs customized to fit your needs.

The SE530 & SE540 sealers can be used stand alone or linked as a multi-head sealing unit.



RAPID SEAL™ SE530 CAT# 428-SE530

- Space saving benchtop sealer
- Fixed sealing head



RAPID SEAL™ SE540 CAT# 428-SE540

- Space saving benchtop sealer
- Fixed and hand-held sealing unit



RAPID SEAL™ SE330 CAT# 428-SE330

- Standard benchtop sealer
- Fixed sealing head



RAPID SEAL™ SE340 CAT# 428-SE340

- Multiple station benchtop sealer
- "Fixed and hand-held sealing unit"



RAPID SEAL™ SE370 CAT# 428-SE370

- Battery powered benchtop sealer
- Fixed and hand-held sealing unit



RAPID SEAL™ SE640 CAT# 428-SE640

- Portable sealing unit
- Hand-held sealing unit



SE800 — Biotechnology & Industrial Applications CAT# 428-SE800

- Ergonomically designed hand-held sealing gun
- Adjustable sealing time handles industrial tubing
- Seals 5 mm-16 mm PVC tubing

CELL PROCESSING



GENESIS™ ES315 ELECTRONIC OPTICAL PLASMA EXPRESSOR

CAT# 422-ES315C

- Audible alarm
- Optical sensor detects red cells
- Automatically clamps; separating plasma from red cells

CORD BLOOD SEALER

CAT# 428-SE470RS

- Specially designed sealing head
- Designed to seal the Pall™ cord blood freezing-bag chamber
- Small footprint / benchtop saving design
- Multiple station use (6' cable)



SEALING JIG

CAT# 428-UCBJIG

- Designed to work with the Pall™ cord blood freezing-bag
- Designed to accurately guide tube sealing head
- Protects freezing bag corners
- Designed to prevent arcing
- Exclusive to GenesisBPS



GENESIS™ CB220

Centrifuge Balance

CAT# 419-CB220

- Electronic differential digital scale
- Displays weight differential of 2 objects
- Balances in seconds
- Simple and easy to operate



To order, call Genesis BPS at

1-866-712-5663